


Isolation of GFP Ensconsin from Bacterial Plasmid and Expression in Cos-7 Cells

Aaron Miller, Shivani Suklikar, and Jeff Bailey

Results



Materials and Methods

Materials. Liquid culture of *E. coli* containing EMTB-3XGFP plasmid. L-broth ampicillin plate. Qiagen plasmid isolation kit from Qiagen (Venlo, Netherlands). 10x reaction buffer C. EcoRI and XhoI restriction enzymes. 1% agarose electrophoresis gel with 1 μ g/mL ethidium bromide. COS-7 cell media consisting of DMEM (Gibco #121600) with 10% fetal bovine serum (Gibco #10438-026) and 1x penicillin-streptomycin-fungizone mixture (Gibco #15240-062), all from Life Technologies (Carlsbad, CA). COS-7 cells (ATCC number: CRL-1651). Trypsin-EDTA. Eugene, Prolong Gold, and DAPI.

Producing colonies of plasmid-containing bacteria. The liquid culture of *E. coli* containing the EMTB-3XGFP plasmid was streaked onto an LB-Amp plate. Plate was incubated at 37°C for

roughly 24 hours.

Isolation of EMTB-3XGFP plasmid DNA. A single colony was used to inoculate 6mL of LB-Amp media and grown for roughly 24 hours. 1.5mL of this liquid culture was micro-centrifuged at 10,000g for 5 minutes; the supernatant was discarded. Another 1.5mL was micro-centrifuged in the same Eppendorf tube. The Qiagen Kit protocol for plasmid purification was followed (QIAprep Miniprep Handbook 12/2006, pg. 22-23).

Analysis of plasmid DNA. The OD₂₆₀ and OD₂₈₀ of the plasmid DNA were measured. Approximately 1 μ g of 150ng/ μ plasmid DNA was put into each of three Eppendorf tubes labelled A, B, and C. Tube A was the negative control, containing only DNA in 1x reaction buffer C; tubes B and C additionally contained restriction enzymes (0.5 μ l EcoRI in the case of tube B and 0.5 μ l of both EcoRI and XhoI in the case of tube C). Four wells of an agarose electrophoresis gel were loaded as follows: 5 μ l ladder, 15 μ l from each of tubes A, B, and C. Gel was run for 47 minutes at 100 volts, and the result was photographed under UV illumination.

Splitting of COS-7 cells. 100mL of COS-7 cell media was prepared. A plate of COS-7 cells was washed with 10mL PBS and trypsinized for 5 minutes. It was verified that the cells were dislodged using an inverted microscope. Trypsin was inhibited with 5mL media, and the cells were collected to a pellet at 1500rpm for 3 minutes; the supernatant was aspirated. Cells were resuspended in media and counted with a hemacytometer. Cells were split to low and high cell count plates.

References

Weimbs, T. (2013). *Cell Biology Laboratory Manual*. MCDB department at University of California, Santa Barbara.

Qiagen Plasmid Purification Handbook. (2006).
<http://www.lifesci.ucsb.edu/mcdb/labs/weimbs/classes/2013-MCDB_103L/additional/QIAprep_Miniprep_Handbook.pdf>.